



# MODULE 1: ULTRASOUND FOUNDATIONS (PHYSICS + INSTRUMENTATION)

Carlos Jimenez PhD (c), DPT, RMSK, CPSS

1

---

---

---

---

---

---

---

---



## Ultrasound Foundations: Physics, Image Optimization, and MSK Relevance

- Why this matters:
- Image clarity
- Artifact recognition
- Safe scanning
- The bridge between **good scanning** and **diagnostic confidence**



2

---

---

---

---

---

---


---

---



## Ultrasound = Mechanical Longitudinal Sound Waves"

- Mechanical, not electromagnetic
- Requires a medium
- Travels as **compression & rarefaction**
- Frequency >20,000 Hz



3

---

---

---

---

---

---

---

---

**Frequency = Resolution vs Penetration**



- ↑ Frequency → ↑ resolution, ↓ depth
- ↓ Frequency → ↓ resolution, ↑ depth

**Clinical examples:**

- 15–18 MHz: tendons, ligaments, nerves
- 6–8 MHz: deep hip, proximal hamstring



4

---

---

---

---

---

---


---

---

**Speed of Sound in Soft Tissue = 1540 m/s**



- Ultrasound machines assume one speed
- Incorrect assumption → artifacts
- Variable tissues:
  - Fat slower
  - Tendon faster
  - Bone fastest (major reflector)



Tissue	Speed of Sound (m/s)
Air	330
Fat	1450
Water	1480
Soft tissue (average)	1540
Liver	1550
Kidney	1560
Blood	1570
Muscle	1580
Bone	4080

5

---

---

---

---

---

---

---


---

**Why Some Structures Look Bright or Dark**



- Impedance = density × speed
- High mismatch → strong reflection (bright)
- Low mismatch → weak reflection (dark)

**Normal Anatomy**



6

---

---

---

---

---

---

---

---



### Specular vs Diffuse Reflection

- Specular (smooth): tendon, nerve, bone cortex
- Diffuse (rough): muscle, fat
- Angle dependency
- Clinical relevance: **anisotropy**




Specular Reflection

One Direction



Diffuse Reflection (Scattering)

Multiple Directions  
Low Amplitude



Diffuse Reflection (Scattering)

Multiple Directions  
Low Amplitude

7

---

---

---

---

---

---

---

---

---

---



### Why Structures Appear Shifted or Duplicated - Refraction

- Occurs at angled boundaries with speed change
- May create "step-off" appearance
- Important in:
- Needle guidance
- MCL, Achilles scans (angled anatomy)

8

---

---

---

---

---

---

---

---


---

---



### Why Images Fade with Depth - Attenuation

- Caused by: absorption, scatter, reflection
- $\uparrow$  Frequency =  $\uparrow$  attenuation
- Compensation tools:
- TGC
- Lower frequency
- Harmonics OFF for deep MSK



Mild      Moderate      Severe

Fatty Liver Examples

9

---

---

---

---

---

---

---

---

---

---

**Lateral, Axial, Elevational Resolution**



How Sharp Your Image Actually Is

- Axial = along beam → best resolution
- Lateral = perpendicular → depends on focusing
- Elevational = slice thickness → bursa vs tendon artifact

10

---

---

---

---

---

---

---

---

**Transducers**

Why Linear Probes Dominate MSK

- Linear probe: high frequency, rectangular field
- Curvilinear: deep hip, glute, proximal hamstring
- Keep the footprint perpendicular

Image	Probe Type	Frequency (MHz)	Use Case
	Linear	10-17 MHz	Superficial structures, tendons, muscles, ligaments (e.g., wrist, elbow, tendon)
	Curvilinear	3-7 MHz	Deeper structures, abdomen, hips
	Compact Linear	15-17 MHz	Small areas, ligaments, dynamic MSK
	Phased Array	2-5 MHz	Cardiac or deep thoracic

11

---

---

---

---

---

---

---

---

**Focusing**

Place the Focal Zone at the Structure of Interest

- Improves lateral resolution
- Narrowest beam = sharpest point
- Multiple focal zones slow frame rate



12

---

---

---

---

---

---

---

---

**The Controls You Actually Need to Master**



- Depth
- Gain
- TGC
- Frequency
- Dynamic Range
- Focus
- Harmonics
- Compound imaging



13

---

---

---

---

---

---


---

---

**Set Depth for the Smallest Field That Fits the Anatomy**



- Too deep = tiny structure
- Too shallow = anatomy cut off
- Ideal: structure centered in screen middle third



14

---

---

---

---

---

---

---

---

**Gain & TGC**

Brightness Should Match Tissue Realism

- Gain = global brightness
- TGC = depth-specific
- Avoid "over-gaining tendons"



15

---

---

---

---

---

---


---

---

**Dynamic Range**

**Contrast Control: How Soft or Sharp**

- High DR = smooth, low contrast
- Low DR = crisp, high contrast
- Useful for: nerves, tendons



16

---

---

---

---

---

---


---

---

**Harmonics & Compound Imaging**

**When Advanced Modes Help — and When They Hurt**

- Harmonics ON: improves clarity in superficial structures
- Harmonics OFF: needed for deeper structures
- Compound imaging: reduces noise, reduces artifacts



17

---

---

---

---

---


---

---

---

**Artifacts Are Not Errors — They Are Clues**

- Anisotropy (most important)
- Shadowing
- Acoustic enhancement
- Reverberation
- Side lobe artifact



18

---

---

---

---

---

---


---

---

 **Anisotropy (Deep Dive)**

**Angle-Dependent Brightness Loss — Not Pathology**

- Occurs with tendons & nerves
- Beam must be perpendicular
- Fix = heel-toe maneuver



19

---

---

---

---

---

---

---

---

 **Ultrasound Safety Principles: ALARA, MI, TI, and Doppler Output**

- Safe → Effective → Evidence-Based MSK Scanning
- Understanding output keeps MSK practice aligned with POCUS & RMSK standards

20

---

---

---

---

---

---

---

---

 **ALARA = As Low As Reasonably Achievable**

- Use lowest output for diagnostic image
- Limit unnecessary dwell time
- Use appropriate presets

21

---

---

---

---

---

---

---

---



### Mechanical Index (MI) — Risk of Cavitation

- MI relates to acoustic pressure + frequency
- High MI = more cavitation potential
- MSK scanning = naturally very low MI

22

---

---

---

---

---

---

---

---



### Thermal Index (TI) — Risk of Tissue Heating

- TI predicts possible temperature rise
- TIS = soft tissue TI
- TIB = bone TI (higher near cortex)
- MSK TI stays very low

23

---

---

---

---

---

---

---

---



### Doppler = Increased Acoustic Power

- Color & spectral Doppler raise MI & TI
- Use Doppler **purposefully and briefly**
- Essential applications in MSK:
  - Neovascularity (Achilles, patellar)
  - Hyperemia
  - Vascular injury assessment

24

---

---

---

---

---

---

---

---



**Simple Rules for Safe & Effective MSK Ultrasound**

- Use appropriate MSK presets (low output, high frequency)
- Keep probe moving unless clinically necessary
- Use Doppler sparingly and purposefully
- Follow ALARA: image quality → minimal power
- Monitor MI/TI but know MSK values stay low

25

---

---

---

---

---

---

---

---



**Module 1 References — Physics, Instrumentation & Safety**

Hedrick WR et al. *Ultrasound Physics and Instrumentation*. Elsevier; 2013.  
 Kremkau FW. *Diagnostic Ultrasound: Principles and Instruments*. Elsevier; 2020.  
 Hoskins PR et al. *Diagnostic Ultrasound: Physics and Equipment*. CRC Press; 2019.  
 Jacobson JA. *Fundamentals of Musculoskeletal Ultrasound*. 3rd ed. Elsevier; 2018.  
 Smith J, Finnoff JT. *PM&R*. 2009;1(1):64–75.  
 AIUM. *Official Statement on Ultrasound Safety*.  
 WFUMB. *Global Guidelines on Ultrasound Safety*.  
 Barnett SB et al. *Safe Use of Ultrasound in Medical Diagnosis*. BIR; 2012.  
 Bassignani MJ. *Ultrasound Clinics*. 2007.  
 • Fessell DP, van Holsbeeck M. *Ultrasound Clinics*. 2007.

26

---

---

---

---

---

---

---

---



**MODULE 2 — Ultrasound Propagation, Tissue Interaction & Artifacts**

Carlos Jimenez PhD (c), DPT, RMSK, CPSS

27

---

---

---

---

---

---

---

---



### Ultrasound Propagation, Tissue Interaction & Artifacts

- Acoustic Impedance
- Reflection / Scattering / Refraction
- Attenuation
- Tissue Echogenicity
- Key MSK Artifacts

28

---

---

---

---

---

---

---

---



### Correct Interpretation Depends on Understanding Wave Behavior

- Image clarity
- Identifying true pathology
- Fixing artifacts
- Accurate tendon, ligament, nerve evaluation

29

---

---

---

---

---

---

---

---



### Acoustic Impedance (Z)

**Why Different Structures Appear Light or Dark**

- $Z = \text{density} \times \text{speed}$
- Bigger mismatch  $\rightarrow$  stronger reflection
- Examples:
  - Bone vs soft tissue = bright cortex + shadow
  - Tendon vs muscle = strong interface echo

Material	Density, $\rho$ (kgm <sup>-3</sup> )	Speed, c (m s <sup>-1</sup> )	Characteristic impedance Z (kgm <sup>-2</sup> s <sup>-1</sup> ) = $\rho c$
Water	1000	1480	1.48
Blood	1060	1570	1.62
Bone	1380-1810	4000	3.75-7.28
Brain	1030	1558	1.55-1.66
Fat	920	1430	1.35
Kidney	1040	1560	1.62
Liver	1060	1570	1.64-1.68
Lung	400	500	0.20
Muscle	1070	1584	1.68-1.74
Tendon	1060	1560	1.65-1.67

30

---

---

---

---

---

---

---

---

 **Reflection (Primary Image Creator)**

**Reflection = The Echoes We Turn Into Images**

- Specular → tendon, bone cortex
- Diffuse → muscle, fat
- Angle-dependent (critical for MSK)

31

---

---

---

---

---

---


---

---

 **Scattering**

**Why Muscle Looks Grainy or Speckled**

- Irregular surfaces
- Multi-direction echoes
- Creates characteristic MSK texture



32

---

---

---

---

---

---


---

---

 **Refraction**

**Beam Bending = Misplaced or Duplicated Structures**

- Occurs when beam crosses angled tissue with different speed
- Can shift structures or create step-off
- Seen around tendons, MCL, Achilles



33

---

---

---

---

---

---

---

---

### Attenuation

**Why Structures Darken With Depth**

- Caused by absorption, scatter, reflection
- Stronger in: muscle → tendon → bone
- Higher frequency = more attenuation

**Fatty Liver Examples**

Images and program courtesy of Dr. Mark Stone, Washington Physicians, PA, Inc.

34

---

---

---

---

---

---

---

---

### Echogenicity Patterns

**How Each MSK Tissue Normally Appears**

- **Bone:** bright cortex, clean shadow
- **Tendon:** fibrillar, bright
- **Muscle:** pennate, hypoechoic with fibroadipose septa
- **Ligament:** linear, compact
- **Nerve:** honeycomb pattern
- **Fat:** hypoechoic with striations

35

---

---

---

---

---

---

---

---

### Bone

**Impedance Wall: Strongest Reflector**

- 99% reflection
- No through-transmission
- Posterior shadow → diagnostic

36

---

---

---

---

---

---


---

---

**Tendon**

Highly Organized, Highly Reflective

- Hyperechoic fibrillar
- Strong angle dependency
- Susceptible to anisotropy



37

---

---

---

---

---

---


---

---

**Muscle**

Pennate Architecture + Scatter + Hypoechoic Tissue

- Hypoechoic muscle fibers
- Hyperechoic septa
- Pennation angle changes with contractio



38

---

---

---

---

---

---

---

---

**Nerves**

Honeycomb in Short Axis, Fascicular in Long Axis

- Fascicles hypoechoic
- Perineurium hyperechoic
- Angle sensitive (mild anisotropy)



39

---

---

---

---

---

---

---

---

 **Ligaments**

**"Compact, Hyperechoic Bands with Linear Fibers"**

**On Slide:**

- Less anisotropic than tendons
- Best seen under tension
- Cortical attachment landmarks help localization



40

---

---

---

---

---

---

---

---

 **Key MSK Artifacts (Overview)**

**Artifacts Are Normal — Learn to Use Them**

- Anisotropy
- Shadowing
- Acoustic enhancement
- Reverberation
- Side lobes
- Edge artifact

41

---

---

---

---

---

---


---

---

 **Anisotropy (Most Important MSK Artifact)**

**Angle-Dependent Brightness Loss**

- Common in tendon, ligament, nerve
- Mimics tears
- Fix: heel-toe to maintain perpendicular angle



42

---

---

---

---

---

---


---

---

**Shadowing**

**When Sound Cannot Penetrate**

- Bone → clean shadow
- Calcification → dirty or clean shadow
- Foreign bodies → variable shadowing



43

---

---

---

---

---

---

---

---

**Acoustic Enhancement**

**Brighter Tissue Deep to Fluid**

- Occurs beneath cysts, bursa fluid, effusion
- Helps distinguish solid vs cystic
- Enhances diagnostic confidence



44

---

---

---

---

---

---

---

---

**Edge Artifact (Refraction at Curved Surfaces)**

**Shadowing at Curved Boundaries**

**On Slide:**

- Seen at tendons, muscles, vessels
- Beam refracts around curved boundary → dropout



45

---

---

---

---

---

---

---


---

**Side Lobe Artifact**



False Echoes from Off-Axis Beams

- Echogenic debris in cysts
- Off-axis beams reflect into main image
- Importance in nerve scanning



46

---

---

---

---

---

---

---

---

**Understanding Tissue Interaction = Accurate Interpretation**



- Predictable tissue appearance
- Artifacts = diagnostic clues
- Angle dependency
- Tissue-specific expectations

47

---

---

---

---

---

---

---

---

**References — Module 2**



- Hedrick WR et al. *Ultrasound Physics and Instrumentation*.
- Kremkau FW. *Diagnostic Ultrasound: Principles & Instruments*.
- Hoskins PR et al. *Diagnostic Ultrasound: Physics & Equipment*.
- Szabo TL. *Ultrasound Imaging: Inside Out*.
- Fessell DP, van Holsbeeck M. *Ultrasound Clinics*.
- Thomas RH. *Semin Musculoskeletal Radiol*.
- Jacobson JA. *Fundamentals of MSK Ultrasound*.
- Bianchi S, Martinoli C. *MSK Ultrasound*.

48

---

---

---

---

---

---

---

---



**MODULE 3 — Pulse Echo  
Instrumentation & Imaging  
Principles**

Carlos Jimenez: PhD (c), DPT, RMSK, CPSS

49

---

---

---

---

---

---

---

---



**Module 3: Pulse Echo  
Instrumentation & Imaging Principles**

- Pulse formation
- Echo reception
- Signal processing
- Scan conversion
- Pre- & post-processing

50

---

---

---

---

---

---

---

---



**SECTION A — PULSE-ECHO  
BASICS**

51

---

---

---

---

---

---

---


---

### What Is Pulse-Echo Imaging?



Transmit → Wait → Receive → Process → Display

- One pulse at a time
- Listening time determines depth
- Echo intensity determines brightness



52

---

---

---

---

---

---

---

---

### Range Equation (Depth Calculation)



- Distance = (Speed × Time) / 2
- Ultrasound assumes 1540 m/s
- Time-of-flight determines depth
- Divided by 2 because sound travels there and back

**Range Equation**

Speed of sound (soft tissue) 1540 m/s    Time (s)

$$D = \frac{C \cdot t}{2}$$

Depth (m)                      Return time (s)

Range equation questions are common on the exam.

53

---

---

---

---

---

---

---


---

### Reflection (Primary Image Creator)



Reflection = The Echoes We Turn Into Images

- Specular → tendon, bone cortex
- Diffuse → muscle, fat
- Angle-dependent (critical for MSK)



54

---

---

---

---

---

---

---

---



## SECTION B — PULSE CHARACTERISTICS

55

---

---

---

---

---

---


---

---



## Pulse Repetition Frequency (PRF)

- How Many Pulses Sent Per Second
- $\uparrow$  PRF =  $\uparrow$  frame rate
- $\downarrow$  PRF = deeper imaging
- High PRF risks range ambiguity



Time (sec)

56

---

---

---

---

---

---


---

---



## Pulse Duration

- Length of Time the Pulse Is 'On'
- Shorter pulses = better axial resolution
- Determined by number of cycles  $\times$  period
- Bandwidth and damping influence this



57

---

---

---

---

---

---

---

---

**Spatial Pulse Length (SPL)**

**Physical Length of the Pulse in Space**

- $SPL = \# \text{ cycles} \times \text{wavelength}$
- Shorter SPL = better axial resolution
- Higher frequency = shorter wavelength

58

---

---

---

---

---

---

---

---

**Duty Factor**

**Percentage of Time Ultrasound Is Transmitting**

- Diagnostic US duty factor <1%
- Echo machine spends 99% of time listening
- Doppler increases duty factor

59

---

---

---

---

---

---

---

---

**SECTION C — THE TRANSMITTER & RECEIVER**

60

---

---

---

---

---

---

---

---

 **Transmitter (Output Section)**

**Controls Pulse Creation**

- Determines:
- Pulse amplitude
- Pulse duration
- PRF

61

---

---

---

---

---

---

---

---

 **Receiver Functions Overview**

**How the Machine Processes Echoes**

- Amplification
- Compensation
- Compression
- Demodulation
- Rejection

62

---

---

---

---

---

---

---

---

 **Amplification (Gain)**

**Boosting All Echo Signals Equally**

- Does **not** increase transmitted power
- Makes entire image brighter
- Too much gain masks pathology




Fig 1 Overall gain set too high      Fig 2 Overall gain set too low

63

---

---

---

---

---

---


---

---

**Compensation (TGC)**

**Corrects for Attenuation With Depth**

- Adjusts brightness at specific depths
- Creates uniform image
- Essential for MSK



64

---

---

---

---

---

---

---

---

**Compression**

**Reduces Wide Range of Echo Intensities**

- Narrows signal dynamic range
- Makes image readable
- User control = dynamic range knob

65

---

---

---

---

---

---


---

---

**Demodulation**

**Prepares Signal for Display**

Two steps:  
1. Rectification  
2. Smoothing (enveloping)



66

---

---

---

---

---

---


---

---

 **Rejection (Threshold)**

**Suppresses Low-Level Noise**

- Eliminates unwanted echoes
- Too much → loss of detail



Amplitudes below dotted line reset to zero

67

---

---

---

---

---

---

---

---

 **SECTION D — SCAN CONVERTER & DIGITAL PROCESSING**

68

---

---

---

---

---

---


---

---

 **Scan Converter**

**Converts Echo Data Into Pixels**

- Receives analog signals
- Stores digital values
- Uses X-Y grid for display



69

---

---

---

---

---

---

---

---



## Preprocessing vs Postprocessing

Machine Decides vs User Decides

- Preprocessing:
  - Dynamic range
  - Persistence
  - Write zoom
- Postprocessing:
  - Map changes
  - Read zoom
  - Gray-scale curve

70

---

---

---

---

---

---

---

---



## Preprocessing

- TGC (Time Gain Compensation)
- Dynamic range (compression ratio)
- Persistence (frame averaging)
- Write zoom
- Edge enhancement (if present)
- Frequency compounding (if selectable)
- Spatial compounding
- RES-level adjustments (line density changes)
- Acoustic power
- Beam steering & focusing adjustments

71

---

---

---

---

---

---

---

---



## Postprocessing

User decides — can be changed anytime after the image is stored or frozen.

These include:

- Gray-scale curve selection
- Map changes / color maps
- Read zoom
- Brightness & contrast display adjustments
- Smoothing filters applied after storage

• These do not change the raw echo data.

72

---

---

---

---

---

---

---

---



## SECTION E — REAL-TIME IMAGING LIMITATIONS

73

---

---

---

---

---

---

---

---



## Frame Rate

**Higher Frame Rate = Better Motion Imaging**

Factors that increase frame rate:

- High PRF
- Narrow sector width
- Low line density
- Shallow depth

74

---

---

---

---

---

---

---

---



## Temporal Resolution

**Ability to Accurately Display Motion**

- Improved by: shallow depth, narrow field, low line density
- Reduced by: high frame averaging, large sector widths

75

---

---

---

---

---

---

---

---



### Things RMSK Will Test

Must-Know Exam Concepts

- Range equation
- Receiver functions
- Pre- vs post-processing
- PRF & depth relationship
- SPL & axial resolution
- Frame rate determinants

76

---

---

---

---

---

---

---

---



### Summary

“Pulse → Process → Display → Interpret”

- Pulse creation
- Echo reception
- Signal processing
- Scan conversion
- Image optimization

77

---

---

---

---

---

---

---

---



### Module 3 References

- Hedrick WR et al. *Ultrasound Physics and Instrumentation*. Elsevier.
- Kremkau FW. *Diagnostic Ultrasound: Principles and Instruments*. Elsevier.
- Hoskins PR et al. *Diagnostic Ultrasound: Physics and Equipment*. CRC Press.
- Szabo TL. *Diagnostic Ultrasound Imaging: Inside Out*. Academic Press.
- AIUM Official Statements on Instrumentation & Image Formation.
- Baker ML et al. *Ultrasound Clinics*. 2014.
- Athanasiou A et al. *Eur J Radiol*. 2012.

78

---

---

---

---

---

---

---

---